

Hemoglobin Meter

Mission Hb[®]

*Easy, Fast and Accurate way
for testing of **hemoglobin***



Accurate and Easy to Use

Fast Results (≤ 15 Seconds) of **Hb** & **HCT**

Capillary & **whole blood** can be used



IVD

Hb METER

Now Testing of
Hemoglobin is easy



- Sample Volume: 10 µl
- Measuring Time: <15 Seconds
- Battery Life: 1000 Tests
- Calibration: Auto Calibration by Code Chip
- Memory: 1000 Tests Results
- Measuring Range: 0-25 gm/dl

Easy 3 Step Procedure



1. Insert Test Strip



2. Apply Sample



3. Read Result



Distributor:

Vortex MedSolutions Pvt Ltd
No.14, Sea Mist, Sector 8, Charkop,
Kandivali West, Mumbai 400067
sales@vortexmeds.com; 022-28675723